

# Forensic Microscopy

## Lab 2: Micrometry

Aritra Bhattacharya  
Graduate teaching assistant  
[ecoevaritra23@tamu.edu](mailto:ecoevaritra23@tamu.edu)

(Grab a Microscope from the black closet)


# Where can you find the slides?



- <https://aritra-evolves.github.io/FIVS215.html>
- Canvas


# Today's Agenda

- ❖ Why Micrometry
- ❖ Taking Measurements
- ❖ Calibrating
- ❖ Kohler Review (if wanted)




# Why micrometry?

- ❖ Proper calibration of the microscopes
- ❖ Accurate measurement of samples
- ❖ Ease in the future



Stage Micrometer





Ocular Micrometer

# Taking Measurements

Two Types of Fibers

Trilobal & Not




# Finding your Dominant Eye

If the object moves, you are looking with your  
NON-DOMINANT eye

If the object does not move, you are looking with  
your DOMINANT





# Steps


1. Focus on the ocular reticle
  - a. Make sure the scale is in your dominant eye, move it if not
2. Swing in the lowest objective
3. Focus on the stage micrometer
4. Align the two micrometers to be parallel
5. Choose two points on the stage micrometer - the farther apart the better
6. Count the number of ocular micrometer marks between the two points
7. Divide using the following equation
8. Repeat with each objective




# Step 1 – Focus the Ocular Micrometer

Make sure the ocular micrometer is in clear focus and level

Turn the eyepiece and adjust the diopters as necessary



**no**



**yes**

# Step 2 – Setting Up the Stage Micrometer


Make sure the lowest objective is swung in (10x)

Place the stage micrometer on the stage and make sure it is in focus

Align the two micrometers so that they are PARALLEL




Rotate the stage dial for alignment



# Step 3 – Count the Ticks

- Select two points on the STAGE micrometer - the farther apart the better
- Count the number of ticks on the OCULAR micrometer between the two points on the STAGE micrometer
- Write down each count



$$1 \text{ Ocular Micrometer Division} = \frac{\# \text{ Stage Micrometer Divisions}}{\# \text{ Ocular Micrometer Divisions}} \times K$$

K = The Value of the STAGE Micrometer Division in MILLIMETERS


K = 0.01mm (for the sake of this example)

1 OMD =  $(60/30) \times 0.01 \text{ mm}$

= 0.02 mm

1 OMD = 20 micrometer (in this example)

1mm = 1000 micrometer




1 DIV = 0.01mm

Value of K



Is it easier to measure  
now?



# Exercise


Set up and calibrate stage and ocular micrometry  
For ALL THE OBJECTIVE LENSES (4X, 10X, 40X, 63X).

Write down your calculated values on the index cards and attach it to your microscope

Keep a copy of the calculated values in your lab notebook.

# See you next week!

Stick around if you want to review Kohler again

A dark blue diagonal gradient bar that starts from the bottom left corner and extends towards the top right corner, covering the lower half of the slide.